



**Jonoke Software Development Inc.**  
8709-102ave  
Edmonton, Alberta  
Canada, T5H 4E5

Ph: (780) 448-3647  
Fax: (780) 448-3741  
www.jonoke.com  
info@jonoke.com

---

## **Installing MediFile\_Server**

**Date of Tech Note: April 4, 2002**  
**Author of Tech Note: Adam Q. Bevan**  
**Program Version: MediFile Version 3.1**  
**Enhancement: N/A**

**Bug Fix: N/A**

### **Description:**

Installing a Windows MediFile\_Server v3

### **Detail**

Note: MediFile Server should only be installed on a Windows NT sp6 or Windows 2000 computer using a manually set TCP/IP network protocol. For further Information on Server Configuration please contact Jonoke at 1-780-448-3647.

1. Put the MediFile v3 Install CD into the computer.
2. Double click on My Computer.
3. Go to the Tools menu and select Folder Options.
4. Click on the View tab
5. Make sure that it is showing hidden files and folders and NOT hiding extensions for known file types.
6. Click on Apply then OK.
7. Double click on the CD Rom drive.
8. Double click on the MediFile v3 folder.
9. Double click on the OEM (dmmmyy) folder.
10. Press control-A then press control-C to copy all files in this folder.
11. Press the Backspace key twice then minimize the OEM window.
12. Double click on My Computer.
13. Double click on the C: drive.
14. Create a new folder and name it Jonoke.
15. Double click on the Jonoke folder.
16. Press control-V and minimize the Jonoke window.
17. Restore the CD Rom volume window.
18. Single click on the 4D Utilities folder and control-click on the Network Components folder.
19. Press control-C and close the CD Rom volume window.
20. Restore the Jonoke window and press control-V.
21. Press control-A and go to the File menu and select Properties.
22. Make sure that all files in the Jonoke folder are NOT read-only.
23. You can create a short cut on the desktop of MediFile\_Server.exe.
24. Copy the Network Components folder and minimize the Jonoke window.
25. Double click on My Computer then double click on the C: drive.
26. Double click on the Windows folder and click on show all files and folders.
27. Create a new folder and name it ACI.

28. Double click on the ACI folder and press control-V.
29. Rename the Network Components folder to Network.
30. Close the ACI window and restore the Jonoke window.
31. Make sure that a printer is setup on the server.
32. Double click on MediFile\_Server.exe.

This should start MediFile Server.

If there are any problems please contact Jonoke.