



Jonoke Software Development Inc.
8709-102ave
Edmonton, Alberta
Canada, T5H 4E5

Ph: (780) 448-3647
Fax: (780) 448-3741
www.jonoke.com
info@jonoke.com

Installing Timbuktu Pro

Date of Tech Note: July 2, 2002

Author of Tech Note: Adam Q. Bevan

Program Version: MediFile Version 3.1.1

Enhancement: N/A

Description:

Installing Timbuktu Pro on Windows 2000

Detail

- Quit MediFile Server
- Insert Timbuktu CD
- Follow installation instructions given by the CD or run the Setup.exe file on the CD and step through each dialog installing Timbuktu in the Programs directory.
- Once finished installing you should restart the server computer.
- Once restarted you should now see a new item in the task bar dock. It should look like two computers with a blue screen(not green). Double click on the new item called Timbuktu Pro.
- Go to the Setup menu and select Security.
- Go to the Guests tab.
- Uncheck all of the check boxes.
- Go to the Registered Users tab.
- Click on the Add Registered User button.
- Type in the Name field "Admin".
- Type in the Password and Confirm field the Admin password.
- Check off Send, Control, Observe, Exchange(read, write, delete).
- Click on the Apply button.
- Type in the Name field "Admin".
- Type in the Password and Confirm field "blue".
- Check off Send, Control, Observe, Exchange(read, write, delete).
- Click on the Apply button.
- Click on the Ok button.
- Go to the Setup menu and select Preferences.
- Go to the General tab.
- Set the computer name as MediFile Server.
- Make sure Start Timbuktu as a service is checked off.
- Go to the Observe/Control tab.
- Make sure Display desktop wallpaper is unchecked.
- Go to the Master Password tab.
- Enter in the Password and Confirm field "4483647".
- Click on the Apply button.
- Click on the Ok button.

Note: Please inform Jonoke of the access to your server and passwords you set for Jonoke. We will change our password the next time we log in to your computer.